

## DAPA Determination Procedure

September 2013

### I. Reference:

Webster, P.M., W.H. Hoover, and T.K. Miller. 1989 Determination of 2, 6 diaminopimelic acid in biological materials using high performance liquid chromatography. Anim. Feed. Sci. Technol., 30:11-20.

### II. Personal Protective Equipment:

1. Lab coat
2. Safety glasses/goggles
3. Designated chemical gloves

### III. Reagents:

1. Mobile Phase A:
  - a)  $\approx$  1600 mL of buffer solution- use this to q.s. to 2 liters
    - Buffer Solution:
      - 1600 mL ddH<sub>2</sub>O
      - 3.2812 g sodium acetate
      - 5.67 g anhydrous sodium phosphate
      - Adjust pH to 6.8 w/ 85% H<sub>2</sub>PO<sub>4</sub>
  - b) 370 mL of methanol
  - c) 30 mL Tetrahydrofuran (THF-**peroxide former**)
2. Mobile Phase B:
  - a) 1300 mL methanol
  - b) 700 mL H<sub>2</sub>O
3. Methanesulfonic acid (4N solution)
  - a) 259 mL methane sulfonic acid
  - b) 2 g 3-(2-aminoethyl) indole (0.2% w/v)
  - c) 100 mL of  $\beta$ -aminoadipic acid stock solution
  - d) Q.s. mixture with ddH<sub>2</sub>O
    - Store at 4°C and flush w/ nitrogen after each use
4. DAPA Stock Solution
  - a) 1 mg DAPA/5mL ddH<sub>2</sub>O
5.  $\beta$ -aminoadipic Acid Stock Solution
  - a) 1 mg  $\beta$ -aminoadipic acid/5mL ddH<sub>2</sub>O

### IV. Working Standard Prep:

- a) Combine 1 mL DAPA stock solution w/ 1 mL  $\beta$ -aminoadipic stock solution
- b) 10 mL of ddH<sub>2</sub>O
- c) 1 mL of working standard solution added to 1 mL of OPA derivatizing solution in small test tube

- d) Vortex contents
- e) Place on HPLC w/ prepped samples

#### V. Sample Prep:

- b) Weigh 0.10-0.20 g of sample into screw cap glass tubes
- c) Add 3 ml 4.0 N methanesulfonic acid containing catalyst and internal standard
- d) Securely cap the tubes (tighten caps and slightly release pressure) and place on hotplate for 22 hours (high setting on hot plate)
- e) Allow samples to cool after 22 hours
- f) Neutralized the solution with 10 N sodium hydroxide-use pH probe and adjust to pH of 7.0-7.5
  - Note if you need to bring the pH down to 7.0 instead of up use 1.0 N sulfuric acid
- g) Sample must then be vortexed to ensure the sample is uniform
- h) 0.2000-0.2500 g of HCl-washed activated charcoal was added to each tube
- i) Contents were mixed and transferred to smaller centrifuge tubes
- j) Centrifuge tubes at 1000 x g for 20 min
- k) Supernatant was passed through a disposable filter assembly (pore size 0.45  $\mu\text{m}$ )
- l) A second clean-up procedure utilizing the Sep-Pak C-18 cartridge (Waters Inc.) was performed on the sample as follows:
  - i) Sep-Pak cartridge was activated with 20 mL methanol, washed with 20 mL of 0.1% TFA in water, followed by 10 mL of 0.1% TFA in a 80:20 mixture of water: methanol
  - ii) 1 mL of sample was mixed with 2 mL of 0.1% TFA in 70:30 water: methanol (sample should be free of particulate matter) and the mixture passed through the Sep-Pak cartridge
  - iii) The first milliliter was discarded and the next 2 mL, which contained all the amino acids, were collected.
- m) 1 mL of OPA derivatizing solution was added to 1.0 mL of sample in a small test tube
- n) Contents were mixed thoroughly and placed on the HPLC