

PEPSIN INSOLUBLE NITROGEN (PIN)
Up-dated September 2013

I. Reference:

Goering, H.K. and P.J. Van Soest. 1970. Forage Fiber Analyses. USDA Handbook No. 379.

II. Personal Protective Equipment:

- A. Lab coat
- B. Safety glasses/goggles
- C. Latex gloves

III. Reagents:

- A. Pepsin - HCl solution
 1. Add 8.5 mL of 12 N HCl carefully to approximately 800 mL distilled water.
 2. Add 10 g pepsin, mix and make up to one liter.

IV. Procedure:

- A. Samples must be freeze-dried to avoid heat damaging the nitrogen in the sample
- B. Weigh 0.5 g dried ground sample into 100 mL in vitro tube
- C. Add 50 ml pepsin-HCl solution
- D. Incubate in water bath at 39 °C for 20 hours
- E. Filter with hot water onto Whatman 541 filter paper. Wash thoroughly.
- F. Run Kjeldahl of filter papers or scrape, mix, and weigh residue for LECO analysis

V. Calculation:

- A. Determine g N in residue:

$$\text{g N} = \%N * \text{residue wt.}$$

- B. Determine PIN as a percent of DM:

$$\%PIN = \frac{\text{g N in residue}}{\text{DM}} \times 100$$

- C. Determine N as a percent of DM:

$$\%N = \frac{\text{g N in residue}}{\text{DM}} \times 100$$

- D. Express PIN a percent of total N

$$\text{PIN, \% total N} = \frac{\text{g N in residue}}{\text{g N}} \times 100$$

$$\frac{\text{g N in residue}}{(\text{Sample wt.})(\text{DM})}$$

$$\frac{\text{g N}}{(\text{Sample wt.})(\text{DM})}$$

$$\frac{\%PIN}{\%N}$$