

# NONPROTEIN NITROGEN OF FEEDSTUFFS

## Updated September 2013

### I. Reference:

Waldo, D.R. and H.K. Goering. 1979. J. Animal. Sci. 49:1560.

### II. Personal Protective Equipment:

- A. Lab coat
- B. Safety glasses/goggles
- C. Latex gloves

### III. Reagents:

- A. 0.15 M Sodium Chloride (NaCl) (8.77 g NaCl/L).
- B. 10% (w/v) TCA ( $C_2HCl_3O_2$ )

### III. Procedure:

- A. Weigh in duplicate or triplicate 2 to 5 g of sample into 250 ml Erlenmeyer flasks, include blanks.
- B. Add 100 to 200 ml .15 M NaCl.
- C. Place on shaking hot plate. Heat to 40°C and shake for 1 hour.
- D. Filter into 250 ml volumetric flasks through S&S Whatman 410 filter paper (unless soluble N is desired) and rinse several times and q.s. with 0.15 M NaCl.
- E. Take 20 ml filtrate and place in 50 ml centrifuge tube containing 20 ml 10% TCA.
- F. Swirl and place in cooler overnight.
- G. Centrifuge at 14,000 x g for 15 minutes.
- H. Digest 10 ml of supernatant via the Kjeldahl procedure.
- I. Rinse digesta into 100 ml volumetric and q.s. with water.
- J. Save a portion for  $NH_3$ -N determination on auto analyzer.

### IV. Calculations:

$$\% \text{ NPN} = \frac{(\text{NH}_3\text{-N, mg/ml} - \text{blank}) (5000)}{(\text{Sample wt}) (\text{DM}) (1000)}$$

5000 = Dilution factor

1000 = Conversion from g to mg

$$\text{NPN, \% of total N} = \frac{\% \text{ NPN}}{\% \text{ Total N}}$$