

## NEUTRAL DETERGENT FIBER

Updated September 2013

### I. References:

Van Soest, P.J., J.B. Robertson and B.A. Lewis. 1991. Methods for dietary fiber, neutral detergent fiber, and non-starch polysaccharides in relation to animal nutrition. *J. Dairy Sci.* 74:3583.

Van Soest, P.J. and W.C. Marcus. 1964. Methods for the determination of cell wall constituents in forages using detergents and the relationship between this fraction and voluntary intake and digestibility. *J. Dairy Sci.* 47:704. (Abstract of 1964 Amer. Dairy Sci. Assoc., Univ. of Ariz. Tucson).

Defines: Hemicellulose, cellulose and lignin

### II. Personal Protective Equipment:

- A. Lab coat
- B. Safety glasses/ goggles
- C. Latex gloves
- D. Thermo pad for removing beakers from hot plate

### III. Procedure:

- A. Weigh out 0.5000-0.5040 g of sample that has been ground through the Wiley mill with a 1-mm screen (**all high starch samples must be ground through the cyclo-tech**). Add 0.5 g of Sodium Sulfite for removal of protein.
- B. Place in a Berzelius 600 ml beaker and add 100 ml of room temperature neutral detergent solution.
- C. Turn the main power lines on and start water through condensers on fiber reflux apparatus.
- D. Place the beakers on the digestion rack (hot plates) and cover with condensers; bring to boil rapidly.
- E. **Record time that boiling begins for each sample.** Allow sample to reflux for one hour.
  - 1. Alpha-amylase may be added in 0.5 mL increments after reflux begins and ten minutes before filtering high starch samples.
- F. While the samples are refluxing:
  - 1. Fill rinse water reservoir and heat to boiling
  - 2. Obtain dry weight on Gooch crucibles or Whatman 541 filter paper (minimum 3 hrs drying time in 100°C oven).
  - 3. Rinse down sides of beaker, with NDF solution, during reflux if sample creeps up side of beaker and is no longer in contact with solution.
  - 4. If the sample boils up condenser, discard sample.
- G. At the end of the reflux period, remove each beaker separately and wipe condenser with moist sponge to remove condensed detergent solution.
- H. Place crucible on manifold and insure a tight seal and strong vacuum. Turn on small amount of suction and add hot water to crucible.
- I. Slowly decant sample solution into crucible or filter. Once bulk of material has passed, open suction valve, wash out beaker with boiling water, and add to crucible. **Make certain that no residue remains in beaker.**
- J. Wash residue with several volumes of boiling water and dry with vacuum.
- K. Dry samples at 60°C for 24 hours or in a 100 C lab oven for a minimum of 6 hours.

L. Transfer to a desiccator, reweigh dry crucible/filter and residue.

III. Calculations:

$$\text{NDF\%} = 100 \times \frac{(\text{Dry Crucible or Filter} + \text{Residue}) - (\text{Dry Crucible or Filter})}{(\text{Sample Wt}) (\text{lab corrected DM})}$$