

ACID DETERGENT FIBER Updated September 2013

I. References:

Van Soest, P.J. 1963. Use of Detergents in the Analysis of Fibrous Feeds. II. A Rapid Method for the Determination of Fiber and Lignin. Journal of A.O.A.C. 46:830.

Defines: Cellulose, Lignin

II. Personal Protective Equipment:

- A. Lab coat
- B. Safety glasses/ goggles
- C. Latex gloves
- D. Thermo pad for removing beakers from hot plate

III. Procedure:

- A. Weigh out 0.5000-0.5040 g of sample that has been ground through the Wiley mill with a 1-mm screen (**high starch samples must be taken through the cyclo-tech**). Add 0.5 g of Sodium Sulfite for removal of protein.
- B. Place in a Berzelius 600 mL beaker and add 100 ml of room temperature acid detergent solution.
- C. Turn the main power lines on and start water through condensers on fiber reflux apparatus.
- D. Place the beakers on the digestion rack (hot plates) and cover with condensers; bring to boil rapidly.
- E. **Record time that boiling begins for each sample.** Allow sample to reflux for one hour.
 - 1. Alpha-amylase may be added in 0.5 mL increments after reflux begins and ten minutes before filtering high starch samples.
- F. While the samples are refluxing:
 - 1. Fill rinse water reservoir and heat to boiling
 - 2. Obtain dry weight on Gooch crucibles or Whatman 541 filter paper (minimum 3 hrs drying time in 100°C oven).
 - 3. Rinse down sides of beaker during reflux if sample creeps up side of beaker and is no longer in contact with solution.
 - 4. If the sample boils up condenser, discard sample.
- G. At the end of the reflux period, remove each beaker separately and wipe condenser with moist sponge to remove condensed detergent solution.
- H. Place crucible on manifold and insure a tight seal and strong vacuum. Turn on small amount of suction and add hot water to crucible.
- I. Slowly decant sample solution into crucible. Once bulk of material has passed, open suction valve, wash out beaker with boiling water, and add to crucible. **Make certain that no residue remains in beaker.**
- J. Wash residue with several volumes of boiling water and dry with vacuum.
- K. Dry samples at 60°C for 24 hours or in a 100 C lab oven for a minimum of 6 hours.
- L. Transfer to a desiccator and reweigh dry crucible/filter and residue.

IV. Calculations:

$$\% \text{ ADF} = 100 \times \frac{(\text{Dry Crucible} + \text{Residue}) - (\text{Dry Crucible})}{(\text{Initial Sample Wt}) (\text{Lab corrected DM})}$$