

IN VITRO DRY MATTER DISAPPEARANCE
Updated September 2013

I. References:

McDougall, E.I. 1948. Biochem. J. 43:99.

Tilley, J.M.A. and R.A. Terry. 1963. J. Br. Grassland Soc. 18:104.

Mertens, D.R. 1993. Rate and extent of digestion. pp 13-51 in Quantitative Aspects of Ruminant Digestion and Metabolism. Forbes and France, Eds. (AB International, Wallingford, United Kingdom).

II. Personal Protective Equipment:

- A. Lab coat
- B. Safety glasses/goggles
- C. Latex gloves

III. Reagents:

- A. 20% (v/v) HCl
- B. 5% Pepsin Solution
 - 1. Add 5 g pepsin to 100 mL volumetric flask
 - 2. q.s. to 100 mL with ddH₂O
- C. McDougall's Buffer:

Compound	g/Liter	g/18 Liter	g/12 Liter	g/9 Liter
Sodium Bicarbonate (NaHCO ₃)	9.8	176.40	117.6	88.2
Sodium Phosphate Dibasic (Na ₂ HPO ₄)	2.44	49.94	33.24	24.97
Potassium Chloride (KCl)	0.57	10.26	6.84	5.13
Sodium Chloride (NaCl)	0.47	8.46	5.64	4.23
Magnesium Sulfate (MgSO ₄ ·7H ₂ O)	0.12	2.16	1.44	1.08
Calcium Chloride (CaCl ₂ ·2H ₂ O)	0.16	2.90	1.92	1.45

***CaCl₂ must be added only after all others are completely in solution.**

- McDougall's Buffer is prepared in 18 liter batches and kept sealed in the walk-in **cooler**. Prior to use a volume is warmed to 39 °C and reduced with a stream of CO₂.
- During the warming and reducing step, urea is added to the McDougall's at the rate of **1.0 gm/liter**.
 *For in vitro ammonia release studies urea is not used.

IV. Sample Preparation:

- A. Weigh approximately 0.5000- 0.5040 g of feed (in triplicate) into 100 ml in vitro tube. Include at least 5 blanks and experiment # 0414 hay standards in triplicate or 3 tubes as

corn cob standards.

V. Inoculum:

- A. Obtain whole rumen contents from 2 fistulated steers, prior to feeding, being fed a 30% concentrate diet.
- B. Squeeze through 4 layers of cheese cloth into a **pre-warmed thermos** (about 2.5 L rumen fluid per square of cheese cloth).
- C. Pour filtered fluid into 500-1000 mL separatory funnels, gas with CO₂, stopper, and place in 39 °C water bath until the particulate matter rises to the top.
- D. Remove lower portion and mix with warmed, reduced McDougall's Buffer (typically 1:1 ratio) with 1 g urea/L of McDougall's buffer.
- E. Maintain buffered inoculum at 39°C and constant CO₂ until used. The time lag before buffer fluid is dispensed should be kept as short as possible, usually no more than 10 minutes.

VI. In vitro Incubation:

- A. When using dry feeds add 1 mL of ddH₂O to each tube. Treat all tubes equally. Randomly distribute tubes to rack in the water bath
- B. Dispense 50 ml of mixture into each tube.
- C. Flush top of tube with CO₂ and place a stopper (No. 6 stopper which has had a 1/16 inch hole drilled through it. Care must be taken to insure that the 1/16 inch hole is not blocked) with a small hole in it in the tube.
- D. Swirl tubes gently to avoid getting sample on the sides of the tube.
- E. Place in water bath at 39 °C for 48 hours, swirling tubes at least twice daily.

VII. Pepsin Incubation:

- A. After 48 hours, remove stoppers and rinse any sample adhering to the stopper into the tube with a **minimum** of distilled water.
- B. Using a repeater add 1 ml of 20% HCl, swirl tubes, add 1 ml of 20% HCl **let rack sit** for a minute add another 1mL of 20% HCl swirl tubes; and finally add 3 ml of 20% HCl swirl tubes (total HCl = 6 ml).
- C. To avoid excess foam, wait until acid has been added to all tubes before swirling.
- D. Add 2 ml 5% pepsin, swirl thoroughly, replace stoppers, and return to water bath.
- E. Label (in pencil) and dry Whatman 541 filter papers (minimum 3 hrs in 100°C oven) and weigh.
- F. After 24-hour pepsin digestion, filter residue with suction through filter paper.
- G. Wash tubes out thoroughly with hot distilled water.
- H. Dry filter paper and residue at 60 °C for 24 hrs or at 100°C for at least 6 hours. Cool in desiccator for no longer than 2-4 minutes and weigh.

VIII. Calculation:

$$\% \text{ IVDMD} = 1 - \frac{[(\text{Residue} + \text{filter paper}) - \text{filter paper}] - \text{blank}}{(\text{Sample wt.})(\text{DM})}$$

$$\text{Blank} = (\text{Blank residue} + \text{filter paper}) - \text{filter paper}$$