

RUMINAL FLUID VFA PROCEDURE

Updated September 2013

I. Reference:

- a. Erwin et al. 1961. J. Dairy Sci. 44:1768 (for luminal fluid preparation).

II. Personal Protective Equipment:

- a. Lab coat
- b. Safety glasses/goggles
- c. Latex gloves

III. Reagents:

- a. 25% (w/v) meta-phosphoric acid
 - i. Determine amount of 25% MPA needed for entire VFA procedure. Will need at least 400mL for the smallest run.
 - ii. Weigh 100g m-phosphoric acid
 - iii. Add 200mL ddH₂O to a 400mL volumetric flask with a stir bar, quantitatively transfer m-phosphoric acid to the flask and dissolve.
 - iv. q.s. to 400mL with ddH₂O
- b. 25 % (w/v) meta-phosphoric acid (H₃PO₄), 25mM 2-EB (internal standard)
 - i. Add about 50mL of 25% m-phosphoric acid to a 100mL volumetric flask.
 - ii. Weigh ca. 0.2904g 2-ethylbutyrate (C₆H₁₂O₂) and quantitatively add to flask. (**Note:** This is a critical step for the samples to be run on the GC. You must get all the 2-EB into the 25% m-phosphoric acid).
 - iii. Note the exact weight of the 2-EB added to calculated molarity.
 - iv. q.s. to 100mL with 25% meta-phosphoric acid (MW of 2-EB is 116.16).
 - v. Place in bottle and shake thoroughly. Store at 4°C.

IV. Sample Preparation:

- a. Freeze and thaw sample (helps to precipitate protein).
- b. Centrifuge at 5,000 x g for 10 minutes.
- c. Pipette 2.0mL supernatant into pre-labeled tube in duplicate.
- d. Add 0.5mL ice cold 25% m-phosphoric acid/2-EB solution to the sample. Vortex tubes after addition of 2-EB solution.
- e. Allow to stand in refrigeration for 30 minutes. After refrigeration centrifuge at 10,000 x g for 15 minutes.
- f. Fill 3mL tuberculin syringe with supernatant without disturbing the pellet.
- g. Filter through filter-tip syringe into GC vial. Fill vial at least ¾ full. (**Note:** use a 25mm GDX Disposable filter, Polethersulfone Filter Media with Polypropylene Housing and a 0.45µm Pore Size).

V. Standard:

a. Stock Standard:

- i. Add the following approximate amounts to a 200mL volumetric flask. Note the exact weights for calculation of actual concentrations.

VFA	Approx. g/200mL	Approx. Stock Conc. mM	MW	Approx. Work Conc. mM
Ac-	0.7500	62.50	60.0	50.0
Pr-	0.2775	18.75	74.0	15.0
Ibut-	0.1101	6.25	88.1	5.00
But-	0.2202	12.5	88.1	10.0
Ival-	0.1021	5.00	102.1	4.00
Val-	0.1021	5.00	102.1	4.00

- ii. q.s. to 200mL with 25% m-phosphoric acid (w/o 2EB).

b. Working Standard:

- i. Pipette 2mL stock std to 8 tubes and add 0.5mL of ice cold 25% m-phosphoric acid and 2-EB solution.
- ii. Vortex tubes and incubate in the refrigerator for 30 minutes.
- iii. Mix standard tubes and transfer to eight GC vials.
- iv. It is important that the standards and the samples have the exact same amount of ISTD (2-EB) in them. Use the same preparation of 25% m-phosphoric acid, 25mM 2-EB for all standards and samples.

VI. Column Specifications and GC Conditions:

a. Column

i. Packing

1. Supelco 11965
 - a. 10% SP-1200
 - b. 1% H₃PO₄
 - c. 80/100 Chromosorb W AW
2. Supelco 12144
 - a. 15% SP-1220
 - b. 1% H₃PO₄
 - c. 100/120 Chromosorb W AW

ii. Measurements

1. Length: 8 feet
 2. I.D.: 2mm
 3. O.D.: ¼ inch
 - iii. Special Instructions
 1. Leave 3 inch space on injector end
 2. Leave 2 inch space on detector end
- b. GC Conditions
 - i. Temperature
 1. Oven: 145°C
 2. Injector: 185°C
 3. Detector (FID): 200°C
 - ii. Flow Rates
 1. Nitrogen (N₂): 20mL/minute
 2. Hydrogen (H₂): 20mL/minute
 3. Compressed Air: 200mL/minute
 - iii. Slope Sens: 0.10
 - iv. Attenuation: 2⁸
- c. Supplies
 - i. Syringe
 1. Hamilton 80358 w/ Hamilton 80458 needle
 2. Hamilton 90393
 - ii. Septum
 1. Thermogreen LB-2 Septa 11.0mm (Supelco 20654)
 - iii. Inlet Liner
 1. PureCol Inlet Liner 20mm (Supelco 20536)
 - iv. Filters
 1. Nitrogen
 - a. Agilent OT3-2
 - b. OMI-2 (Supelco 23906)
 2. Hydrogen
 - a. Agilent OT3-2
 - b. OMI-2 (Supelco 23906)
 3. Air
 - a. Supelco 20618 w/ Supelco 20298 replacement mesh