

XYLOSE DETERMINATION
Updated September 2013

I. Reference:

A. P. Trinder. 1975. Analyst. 100:12.

II. Reagents:

A. 6.3682 N HCl

1. 52.63% HCl (v/v)

B. Color Reagent

1. Add 0.5 g Phloroglucinol ($C_6H_6O_3$)

2. q.s. to 100 mL with glacial acetic acid ($C_2H_4O_2$)

C. 0.2% Benzoic Acid ($C_7H_6O_2$)

1. Add 2.0 g Benzoic Acid to 1L volumetric

2. q.s. to 1L with dH_2O

3. Heat and stir on stir plate

NOTE: Handle chemicals in designated fume hood.

III. Personal Protective Equipment:

A. Lab coat

B. Safety glasses/goggles

C. Latex gloves

IV. Procedure:

A. Add 0.95 ml 6.3682 N HCl to test tube.

B. Add 50 μ l deproteinized blood (H_2O or blood w/o xylose ($C_5H_{10}O_5$))

C. Vortex thoroughly

D. Add 5 ml color reagent and mix

E. Place in boiling water bath for exactly 5 min

F. Cool rapidly

G. Read and record absorbance at 554 nm

V. Standards:

A. Stock standard

1. Add .5 g Xylose to a 100mL volumetric

2. q.s. to 100 mL with .2% Benzoic Acid

Mg/dL

5

10

15

20

25

30

35

40

ML stock/100 mL .2% Benzoic Acid

1

2

3

4

5

6

7

8
