

COUNTING RUMINAL PROTOZOA

Updated September 2013

I. References:

- Boyne, A. W., J. M. Eadie and K. Raitt. 1957. J Gen Microbiol. 17:414.
Dehority, B. A. 1984. Appl Environ Microbiol. 48:182.

II. Reagents:

A. Phosphate buffer:

1. Sodium Chloride (NaCl) 0.50%
2. Anhydrous sodium acetate ($C_2NaO_2H_3$) 0.13%
3. Potassium phosphate monobasic (KH_2PO_4) 0.03%
4. Potassium phosphate dibasic (K_2HPO_4) 0.10%
5. Magnesium Sulfate ($MgSO_4 \cdot 7H_2O$) 0.01%

a. Adjust pH between 6.9 and 7.2

B. 10% formalin solution

C. Glycerol-Buffer-methyl green mixture:

1. Phosphate buffer 70%
 2. Glycerol 30%
- a. Place on a stir plate and mix. While mixing, add methyl green dye powder sparingly until desired color is reached.
- b. If glycerol solution is too light expect poor staining.

III. Ruminal Sampling:

A. Ruminal protozoa can be counted per gram of ruminal contents or per ml of ruminal fluid.

1. For counts per gram of ruminal fluid.
 - a. Ruminal contents are hand mixed as thoroughly as possible and a 50 ml plastic beaker is used to collect the sample.
 - b. If hand mixing is not possible, composite samples of contents from various locations (reticulum, cranial sac, dorsal sac, and ventral sac).
2. For counts per ml of ruminal fluid.
 - a. Use the suction strainer device to sample ruminal fluid in a 50 or 100 ml plastic tube.

IV. Dilution of Ruminal Contents:

A. **For counts per gram of ruminal fluid.**

1. Pipette approximately 10 ml of ruminal contents with a wide mouth pipette into a tared beaker and note the weight of the sample.
2. To the beaker, pipette 10% formalin solution with a pasteur pipette to equal the weight of ruminal contents (1:1 dilution).
3. A 2.0 ml aliquot of the well mixed, formalinized sample is pipetted with a wide orifice pipette into a test tube containing 18 ml glycerol-buffer-methyl green mixture (1:20 final dilution).
4. The mixture is allowed to stand for at least 4 h before counting, overnight results in better staining (After protozoa are fixed, stain within one week).
5. Make preparations in duplicate.

B. **For counts per ml of ruminal fluid.**

1. Add 5 ml of ruminal fluid to 5 ml of formalin solution (dilution 1:1) and mix.
2. Then place a 2.0 ml aliquot of the well mixed, formalinized sample with a wide orifice pipette into a test tube containing 18 ml glycerol-buffer-methyl green mixture (1:20 final dilution).
3. The mixture is allowed to stand for at least 4 h before counting, overnight results in better staining (After protozoa are fixed, stain within one week).
4. Make tubes in duplicated.

V. Counting the protozoa:

- A. Mix the tube well immediately before counting.
- B. Fill the counting chamber (Sedgewick-Rafter Cell) and place a coverslip on the chamber carefully, without causing formation of air bubbles.
- C. The slide is set for 10 minutes to allow settling of protozoa.
- D. If dilution is too concentrated, further dilution with phosphate buffer prior to counting works well (7-8 protozoa per field is desirable).
- E. Use the eyepiece to count the number of organisms.
- F. Count at least 25 microscopic fields using a 10X objective.
 1. It is better to count 50 fields through all the chambers.

VI. Constants and Calculations (for 50 fields):

A. **Constants:**

1. Eye piece standard square units (one field): 0.469225 mm^2
2. Depth of Sedgewick-Rafter cell: 1.00 mm
3. Standard volumetric units: 0.469225 mm^3
4. Volume of 50 field: 23.461250 mm^3
5. Organisms per count ($1000 \text{ mm}^3 / 23.46125 \text{ mm}^3$): 46.6235

VII. Calculations:

Number of organisms in 50 fields x dilution x organisms per count = Organisms/ml
or g

Example: $184 (50 \text{ fields}) \times 20 \times 46.6235 = 156,854/\text{ml or g}$

Note: For identification of species use an atlas of microbiology. Practice the procedure before you experiment. Be consistent on time of sampling relative to feeding time.